

# The High-pO<sub>2</sub> Method to Prevent U-Loss During Single-Aliquot Hematite (U-Th)/He Measurement: Development, Implementation, and Automation

Florian Hofmann<sup>1</sup>

<sup>1</sup>University of Southern California

November 23, 2022

## Abstract

Single-aliquot (U-Th)/He dating of iron-oxides requires less mass and has a higher spatial resolution than the two-aliquot approach, and is, therefore, a more reliable tool for quantifying the timescales of weathering processes, fault activity, and the development of soils and surfaces. Highly helium-retentive hematite samples must be heated to 1000-1100 °C to completely degas He, but we show that U is lost progressively during laser-heating starting at ~900 °C, with major U-loss at ~980 °C and complete U-loss at 1050-1100 °C. This partial or complete loss of U leads to incorrect (U-Th)/He ages that appear older than the true ages. We performed a series of heating tests on hematite and goethite aliquots with independently determined (U-Th)/He ages of 10-1761 Ma in which the helium release was determined by isotope-dilution mass spectrometry and phase change was monitored by infrared spectroscopy. As hematite is heated, reduction of Fe<sup>3+</sup> causes a phase change to magnetite and then to elemental Fe. We correlate the onset of U-loss to the phase change from hematite to magnetite. By raising the temperature of this phase transition using a high oxygen partial pressure (pO<sub>2</sub>) in the sample chamber during laser-heating, we show that the onset temperature of U-loss is equally raised. Samples can therefore be safely degassed at higher temperatures without any detectable U-loss. In our implementation of this technique, the O<sub>2</sub> in the sample chamber is being withdrawn from a tank using a pipette and it is being released before and captured after the degassing process on activated charcoal in a cold finger with a movable LN<sub>2</sub> Dewar flask. We describe the automation of this process for routine degassing of hematite samples. We show that an average age calculated on a reference hematite sample from replicate aliquots (n=12), which were analyzed using this procedure, has a relative uncertainty of 2% (1σ), and is within uncertainty of the previously measured two-aliquot age. We suggest this high-pO<sub>2</sub> degassing procedure as a way to precisely and reproducibly determine single-aliquot hematite and goethite (U-Th)/He ages.

The high-pO<sub>2</sub> method to prevent U-loss during single-  
aliquot hematite (U-Th)/He measurement:  
Development, implementation, and automation

Florian Hofmann  
University of Southern California


Jonathan Treffkorn  
Ken Farley  
California Institute of Technology

Thermo2021  
September 13, 2021

# Utility of hematite (U-Th)/He


## Geochronology

Pidgeon et al. (2004) – formation age of ferruginous nodules from lateritic duricrusts




## Thermochronology


Farley and Flowers (2013) – temperature history modeling of polycrystalline aggregates



Cooper et al. (2016) – downward migration of water table resulting from incision




Ault et al. (2015) – resetting of hematite crystallites due to shear-heating




# Loss of U during laser-heating


## Ideal case



## Problematic case




# Loss of U during laser-heating





# Natural variability in hematite and goethite samples




## Hematite




## Goethite



# Heating tests of multiple samples




# No difference between Pt and Nb tubes



# What happens during laser-heating?




Reduction of hematite:  
 $3\text{Fe}_2\text{O}_3 \rightarrow 2\text{Fe}_3\text{O}_4 + \text{O}_2$



Phase diagram after  
Ketteler et al. (2001)


# Monitoring phase changes during laser-heating with FTIR

Calibration of FTIR spectra  
with known mixtures of hematite and magnetite




FTIR measurements of heated hematite samples


vacuum, 10 min




vacuum, 10 min



# Hematite to magnetite transition detected by pressure increase




Delaying phase change (and U-loss?) to higher temperatures




Phase diagram after Ketteler et al. (2001)

# Data from laser-heating experiments vs. prediction







- Laser-heating of Pt packets
- FTIR on powdered sample



Hofmann et al. (2020)

Phase diagram after  
Ketteler et al. (2001)

# Effect of increasing pO<sub>2</sub> on U-loss temperature




# Correlation of U-loss and phase change


Delay of phase change/U-loss:

~980 °C → 1200-1300 °C



with ~100 Torr (130 mbar) O<sub>2</sub>




Hofmann et al. (2020)




# Modifications to the quad line to automate this method




# High-pO<sub>2</sub> degassing procedure



# High-pO<sub>2</sub> degassing procedure




# High-pO<sub>2</sub> degassing procedure





# Isothermal holding at 1000 °C in vacuum and O<sub>2</sub>






- Loss of U in vacuum increases with holding time
- No detectable loss of U with high-pO<sub>2</sub> at 4 h of holding



# Analysis of highly-retentive hematite sample using high-pO<sub>2</sub> method




# Recent paper on weathering and groundwater in the Andes using this method



# Conclusions

- Hematite and goethite show detectable U-loss at  $\sim 980$  °C, massive loss at 1050-1100 °C  
→ Pyrometric feedback during laser-heating is important!
- Th/U or Sm/U can indicate major loss of U, but intra-sample variability (usually 10-20%) can make it hard to detect small but significant loss
- Some hematites must be heated to 1000-1100 °C to be completely degassed
- U-loss correlates with phase change from hematite to magnetite
- This phase change and U-loss can be delayed to higher temperatures with an increased oxygen partial pressure ( $pO_2$ ) during laser-heating
- This procedure can be automated for routine hematite (U-Th)/He dating



Chemical Geology 532 (2020) 119350

Contents lists available at ScienceDirect



Chemical Geology

journal homepage: [www.elsevier.com/locate/chemgeo](http://www.elsevier.com/locate/chemgeo)



U-loss associated with laser-heating of hematite and goethite in vacuum during (U-Th)/He dating and prevention using high O<sub>2</sub> partial pressure

Florian Hofmann\*, Jonathan Treffkorn, Kenneth A. Farley

California Institute of Technology, 1200 E California Blvd, Pasadena, CA 91125, United States

