


Optimizing $^{40}\text{Ar}/^{39}\text{Ar}$ analyses using an Isotopx NGX-600 mass spectrometer

INSTRUMENT FEATURES

- Customizable collector block which can consist of a combination of Faraday cup and ion counting electron multiplier detectors.
- ATONA technology performs similarly to 10^{13} - $10^{14} \Omega$ RTIA amplifiers (Cox et al., 2020).

PROJECT GOALS

- Quantify sensitivity of the ATONA backed Faradays.
- Determine optimal measurement and integration durations for various signal sizes.
- Develop and apply a mass discrimination correction.



L5	L4*	L3	L2	L1	Ax	H1	H2	H3	H4
^{36}Ar		^{37}Ar		^{38}Ar		^{39}Ar		^{40}Ar	
	^{36}Ar		^{37}Ar		^{38}Ar		^{39}Ar		^{40}Ar

*IC Multiplier in L4 position. Faradays in all other positions

Sensitivity


Instrument	Moles/fA $^{40}\text{Ar}^*$
ARGUS V	7.00E-18
ARGUS VI	4.00E-18
Nu Noblesse	6.45E-17
MAP215-50	1.13E-17
Helix MC Plus	3.50E-17
NGX-600	$1.21 \pm 0.13 \text{ E-16}$

Differences in sensitivity are dependent on resolution and trap current.

A new cathode source installed on Dec. 3rd can achieve a trap current of up to 1000 μA , which we expect will boost the sensitivity.

Optimization

Longer integration times are ideal (e.g. 100 s). As a result, longer durations are also necessary (e.g. > 800 s) for signal sizes < 0.1V.



Mass Discrimination Correction

Three corrections must be made:

- 1) ATONA gain quantifiable via Isotopx internal calibration.
- 2) Inter-Faraday bias shown to be negligible via peak hopping routine.
- 3) All-Faraday configuration allows for exponential correction using air.

$$\frac{^{40}\text{Ar}}{^{36}\text{Ar}}_{\text{corrected}} = \left(\frac{^{40}\text{Ar}}{^{36}\text{Ar}} \right)_m \times \left(\frac{40}{36} \right)^\beta$$

$$\text{where } \beta = \left(\frac{\ln \left(\frac{(^{40}\text{Ar}/^{36}\text{Ar})_a}{(^{40}\text{Ar}/^{36}\text{Ar})_m} \right)}{\ln \left(\frac{40}{36} \right)} \right)$$

Key Takeaways

- ***Sensitivity*** on cathode source is expected to improve. Advantages of the NGX-600 include the high resolution, large dynamic range, and all-Faraday with ATONA configuration.
- ***Optimizing*** measurement conditions for signals of various sizes and carefully characterizing blanks and baselines is fundamentally important to making best practice measurements.
- ***Mass bias corrections*** are simplified with measurements on a single collector type.

Next Steps

- Define threshold for ^{36}Ar that can be accurately measured on an ATONA backed Faraday
- Measure a range of previously characterized geologic samples (single crystal fusions and groundmass, ranging from Permian to Pleistocene in age).
- Compare current source with Isotopx new cathode source.

Thank you!

Happy to take questions!

Emily Mixon: eemixon@wisc.edu

Acknowledgements: This material is based upon work supported by the National Science Foundation Graduate Research Fellowship Program awarded to E.E.M under Grant No. DGE-1747503. Special thanks to Isotopx for installation and support, to Peter Sobol at UW Madison for assisting with the extraction line build, to Allen Schaeen for formative discussions, and to Jake Ross for endless support with Pychron!

